



# ProTrack™ Pigtail Wire



**FLEXIBLE** Tip  
**STIFF** body



# ProTrack™ Pigtail Wire

## SAFETY — Flexible Spiral Tip


- ▶ Reduce risk of perforation
- ▶ Prevent loss of left atrium access

## SUPPORT — Stiff Tapered Body

- ▶ Confidently push transseptal sheath across difficult septum
- ▶ Facilitate exchange for transseptal sheath with large OD

## Specifications

Model	Body Outer Diameter	Overall Length	Functional Exchange Length
PTW-25-175	0.025 in	175 cm	152 cm
PTW-25-230	0.025 in	230 cm	207 cm



© Copyright Baylis Medical Company Inc., 2013. Baylis Medical Company Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered. ProTrack™ and the Baylis Medical logo are trademarks or registered trademarks of Baylis Medical Company Inc. in the United States of America and/or other countries.

Before use, consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

CAUTION: Federal Law (USA) restricts the use of this device to or by the order of a physician.

Patents Pending and/or issued  
CAR1046

**Baylis**  
MEDICAL

**Baylis Medical Company Inc.**  
5959 Trans-Canada Highway  
Montreal, QC Canada H4T 1A1

Tel. : (514) 488-9801 / Fax : (514) 488-7209  
[www.baylismedical.com](http://www.baylismedical.com) / [info@baylismedical.com](mailto:info@baylismedical.com)